


G6PD EQA MIS System Howto

1. Connect to ***https:// g6pd.qap.tw /MIS/***
2. Use default account : " _____"/password : " _____"to login (You could click "***Change Password***" to change default password. If thus, you should login by changed password)

*** If there is a problem with the website’s security certificate, please click continues.**



Receipt reporting

Clicking "***Sample Receipt***" item on the left of page, you will enter receipt-reporting page. Fill all necessary fields and click submit (If there is any blank field , it will display a warning message). The receiving information will be shown on next screen and wait for your confirmation. You could go back to modify input by clicking [click here to previous page]. After checking your input, click [Confirmed]. If showing “Report sample receiving time successfully”, you have finished this process.

External Quality Assurance for Neonatal G6PD Screening MIS System

G6PD QC Sample Receipt

Your are login as
A03

- [Sample Receipt](#)
- [Input Result](#)
- [View Result](#)
- [Change Password](#)
- [<Online Note>](#)
- [<Back to home>](#)

Screening unit:	A03 xxxxxx
Survey No.:	NS2016-09
Sample set on:	2016-08-22
Sample regular receiving date:	2016-08-23
Report deadline:	2016-08-25
Sample receiving on:	2016 - 08 - 22 PM 4 ←
Also receive next EQA survey QC samples:	NS2016-11 <input checked="" type="checkbox"/> Yes ←
Temperature indicator:	highest temperature which turn black None ←
Humidity indicator:	highest humidity which turn pink 10% ←

External Quality Assurance for Neonatal G6PD Screening MIS System

G6PD QC Sample Receipt

Your are login as
A03

- [Sample Receipt](#)
- [Input Result](#)
- [View Result](#)
- [Change Password](#)
- [<Online Note>](#)
- [<Back to home>](#)

Screening unit:	A03 xxxxxx
Survey No.:	NS2016-09
Sample set on:	2016-08-22
Sample regular receiving date:	2016-08-23
Sample receiving on:	2016-08-22 PM 4
Also receive next EQA survey QC samples:	NS2016-11 Yes
Temperature indicator:	None
Humidity indicator:	10%

Please confirm your input result, you could [click here to previous page](#)

External Quality Assurance for Neonatal G6PD Screening MIS System

G6PD QC Sample Receipt

Your are login as
A03

- [Sample Receipt](#)
- [Input Result](#)
- [View Result](#)
- [Change Password](#)
- [<Online Note>](#)
- [<Back to home>](#)

Screening unit:	A03 xxxxxx
Survey No.:	NS2016-09
Sample set on:	2016-08-22
Sample regular receiving date:	2016-08-23
Sample receiving on:	2016-08-22 PM 4
Also receive next EQA survey QC samples:	NS2016-11 Yes
Temperature indicator:	None
Humidity indicator:	10%
Reporting date:	2016-08-22

Report sample receiving time successfully!
(You don't need to fax this page!)

Result reporting

Clicking "**Input Result**" item on the left of page, you will enter result reporting page (Notice: You could run this process only after you reported sample receiving time). Fill all necessary fields and click [Submit] (Notice: if you select " laboratory prepared " for " Reagent ", the fields of " Assay kit ", " Catalog No ", " Reagent Lot No " could be blank). The input of analysis result will be shown on next screen and wait for your confirmation. You could go back to modify input by clicking [click here to previous page]. After checking your input, click [Confirmed]. If showing "Result input successfully", you have finished this process. You could click "**View Result**" item to view your input of analysis result.

PS. If there is any blank field, it will display a warning message and you should go back to fill the blank field in previous page.

* Any suggestion or comment are welcome, please email to

g6pd@g6pd.tw